

Table of Contents

1 Welcome to eConEXG documentation page.

- [1.1 Installation](#)
- [1.2 Quick Start](#)
- [1.3 Platform Support](#)
- [1.4 More](#)

I Devices

2 iRecorder

- [2.1](#)
 - [2.1.1](#)
 - [2.1.2](#)
 - [2.1.3](#)
 - [2.1.4](#)
 - [2.1.5](#)
 - [2.1.6](#)
 - [2.1.7](#)
 - [2.1.8](#)
 - [2.1.9](#)
 - [2.1.10](#)
 - [2.1.11](#)
 - [2.1.12](#)
 - [2.1.13](#)
 - [2.1.14](#)
 - [2.1.15](#)
 - [2.1.16](#)
 - [2.1.17](#)
 - [2.1.18](#)
 - [2.1.19](#)
 - [2.1.20](#)
 - [2.1.21](#)
 - [2.1.22](#)

3 iFocus

- [3.1](#)
 - [3.1.1](#)

- 3.1.2
- 3.1.3
- 3.1.4
- 3.1.5
- 3.1.6
- 3.1.7
- 3.1.8
- 3.1.9
- 3.1.10
- 3.1.11
- 3.1.12
- 3.1.13
- 3.1.14
- 3.1.15

I.I Trigger Box

4 Wireless

- 4.1
- 4.2

5 Wired

- 5.1
- 5.2

6 Light Stimulator

- 6.1
- 6.2

7 Changelog

- 7.1 0.1.18
- 7.2 0.1.15
- 7.3 0.1.14
- 7.4 0.1.13
- 7.5 0.1.12
- 7.6 0.1.11
- 7.7 0.1.10
- 7.8 0.1.09

1 Welcome to eConEXG documentation page.

eConEXG is a Python SDK for hardwares made by [Niantong Intelligence](#).

1.1 Installation

```
pip install econexg
```

Optional dependencies can be installed via `pip install econexg[option1,option2,...]`

Available options:

bdf: save to BDF file functionality support.
lsl: LSL stream functionality support.
wifi: iRecorder W32 wifi model support, msvc build tools required.

1.2 Quick Start

Example code can be found in [examples](#).

1.3 Platform Support

Hardware	Windows	macOS	Linux
iRecorder USB	✓	✓	✓
iRecorder W8 16	✓	✗	✗
iRecorder W32	✓	✓	✓
iFocus	✓	✓	✓
TriggerBox	✓	✓	✓

1.4 More

- [Discussion Forum](#), to discuss the project on GitHub
- [Issue Tracker](#), if you run into bugs or have suggestions

I. Devices

2 iRecorder

2.1 iRecorder

Bases: Thread

Source code in `src\ConEXG\iRecorder\device.py`

```

13  class iRecorder(Thread):
14      class Dev(Enum):
15          SIGNAL = 10 # signal transmission mode
16          SIGNAL_START = 11
17          IMPEDANCE = 20 # impedance transmission mode
18          IMPEDANCE_START = 21
19          IDLE = 30 # idle mode
20          IDLE_START = 31
21          TERMINATE = 40 # Init state
22          TERMINATE_START = 41
23
24      def __init__(
25          self,
26          dev_type: str
27      ):
28          """
29          Args:
30              dev_type: iRecorder device type. available options: Literal["W8", "USB8", "W16",
31 "USB16", "W32", "USB32"]
32          Raises:
33              Exception: if device type not supported.
34              Exception: if adapter not available.
35          """
36          if dev_type not in {"W8", "USB8", "W16", "USB16", "W32", "USB32"}:
37              raise ValueError("Unsupported device type.")
38          super().__init__(daemon=True, name=f"iRecorder {dev_type}")
39          self.handler = None
40          self.__info_q = Queue(128)
41          self.__with_q = True
42          self.__error_message = "Device not connected, please connect first."
43          self.__save_data = Queue()
44          self.__status = iRecorder.Dev.TERMINATE
45          self.__lsl_flag = False
46          self.__bdf_flag = False
47          self.__dev_args = {"type": dev_type}
48          self.__dev_args.update({"channel": self.__get_chs()})
49
50          self.__parser = Parser(self.__dev_args["channel"])
51          self.__interface = get_interface(dev_type)(self.__info_q)
52          self.__dev_sock = get_sock(dev_type)
53          self.__dev_args.update({"AdapterInfo": self.__interface.interface})
54
55          self._bdf_file = None
56          self.dev = None
57
58          self.set_frequency()
59          self.update_channels()
60
61      def find_devs(self, duration: Optional[int] = None) -> Optional[list]:
62          """
63          Search for available devices, can only be called once per instance.
64
65          Args:
66              duration: Search interval in seconds, blocks for about `duration` seconds and
67          return found devices,
68                  if set to `None`, return `None` immediately, devices can later be acquired by
69          calling `get_devs()` in a loop.
70
71          Returns:
72              Available devices.
73
74          Raises:

```

```

75             Exception: If search thread already running or iRecorder already connected.
76         """
77         if self.is_alive():
78             raise Exception("iRecorder already connected.")
79         if self.__interface.is_alive():
80             raise Exception("Search thread already running.")
81         self.__interface.start()
82         if duration is None:
83             return
84         start = time.time()
85         while time.time() - start < duration:
86             time.sleep(0.5)
87         self.__finish_search()
88         return self.get_devs()
89
90     def get_devs(self, verbose: bool = False) -> list:
91         """
92             Get available devices. This can be called after `find_devs(duration = None)` in a
93             loop,
94             each call will *only* return newly found devices.
95
96             Args:
97                 verbose: if True, return all available devices information, otherwise only return
98                 names for connection,
99                 if you don't know what this parameter does, just leave it at its default
100                value.
101
102             Returns:
103                 Newly found devices.
104
105             Raises:
106                 Exception: adapter not found or not enabled etc.
107             """
108             ret = []
109             time.sleep(0.1)
110             while not self.__info_q.empty():
111                 info = self.__info_q.get()
112                 if isinstance(info, list):
113                     ret.append(info if verbose else info[-1])
114                 elif isinstance(info, bool):
115                     if verbose:
116                         ret.append(info)
117                     elif isinstance(info, str):
118                         raise Exception(info)
119             return ret
120
121     def get_dev_status(self) -> str:
122         """
123             Get current device status.
124
125             Returns:
126                 "SIGNAL": data acquisition mode
127                 "IMPEDANCE": impedance acquisition mode
128                 "IDLE": idle mode
129                 "TERMINATE": device not connected or connection closed.
130             """
131             return self.__status.name
132
133     def get_dev_info(self) -> dict:
134         """
135             Get current device information, including device name, hardware channel number,
136             acquired channels, sample frequency, etc.
137

```

```

138     Returns:
139         A dictionary containing device information, which includes:
140             `type`: hardware type;
141             `channel`: hardware channel number;
142             `AdapterInfo`: adapter used for connection;
143             `fs`: sample frequency in Hz;
144             `ch_info`: channel dictionary, including channel index and name, can be
145 altered by `update_channels()`.

146     """
147     return deepcopy(self.__dev_args)

148

149     @staticmethod
150     def get_available_frequency(dev_type: str) -> list:
151         """Get available sample frequencies of different device types.

152

153     Returns:
154         Available sample frequencies in Hz.
155     """
156     if "USB" in dev_type:
157         return [500, 1000, 2000]
158     return [500]

159

160     def set_frequency(self, fs: Optional[int] = None):
161         """Update device sample frequency, this method should be invoked before
162 `connect_device`.

163

164     Args:
165         fs: sample frequency in Hz, if `fs` is set to `None` or not in
166 `get_available_frequency()`,
167             it will fall back to the lowest available frequency.

168

169     Raises:
170         Exception: Device is already connected.

171

172     New in:
173         - now you can set the sample frequency after device connection.

174     """
175     if self.__status not in [iRecorder.Dev.IDLE, iRecorder.Dev.TERMINATE]:
176         warn = "Device acquisition in progress, please `stop_acquisition()` first."
177         raise Exception(warn)
178     available = self.get_available_frequency(self.__dev_args["type"])
179     default = available[0]
180     if fs is None:
181         fs = default
182     if fs not in available:
183         print(f"Invalid sample frequency, fallback to {default}Hz")
184         fs = default
185     self.__dev_args.update({"fs": fs})
186     self.__parser._update_fs(fs)
187     if self.dev is not None:
188         self.dev.set_fs(fs)

189

190     def connect_device(self, addr: str) -> None:
191         """
192             Connect to device by address, block until connection is established or failed.

193

194     Args:
195         addr: device address.

196

197     Raises:
198         Exception: if device already connected or connection establishment failed.

199     """
200     if self.is_alive():

```

```

201         raise Exception("iRecorder already connected")
202     try:
203         ret = self.__interface.connect(addr)
204         self.__dev_args.update({"name": addr, "sock": ret})
205         self.__dev_args.update({_length": self.__parser.packet_len})
206         self.dev = self.__dev_sock(self.__dev_args)
207         self.__parser.batt_val = self.dev.send_heartbeat()
208         self.__error_message = None
209         self.__status = iRecorder.Dev.IDLE_START
210         self.start()
211     except Exception as e:
212         self.__error_message = "Device connection failed."
213         self.__finish_search()
214         raise e
215
216     def update_channels(self, channels: Optional[dict] = None):
217         """
218             Update channels to acquire, invoke this method when device is not acquiring data or
219             impedance.
220
221             Args:
222                 channels: channel number and name mapping, e.g. `{{0: "FPz", 1: "Oz", 2: "CPz"}`,
223                             if `None` is given, reset to all available channels with default names.
224
225             Raises:
226                 Exception: if data/impedance acquisition in progress.
227         """
228         if self.__status not in [iRecorder.Dev.IDLE, iRecorder.Dev.TERMINATE]:
229             warn = "Device acquisition in progress, please stop_acquisition() first."
230             raise Exception(warn)
231         if channels is None:
232             from .default_config import getChannels
233
234             channels = getChannels(self.__dev_args["channel"])
235             self.__dev_args.update({"ch_info": channels})
236             ch_idx = [i for i in channels.keys()]
237             self.__parser._update_chs(ch_idx)
238
239     def start_acquisition_data(self, with_q: bool = True):
240         """
241             Send data acquisition command to device, block until data acquisition started or
242             failed.
243
244             Args:
245                 with_q: if True, signal data will be stored in a queue and **should** be acquired
246                         by calling `get_data()` in a loop in case data queue is full.
247                         if False, new data will not be directly available and can only be acquired
248                         through `open_lsl_stream` and `save_bdf_file`.
249
250             Raises:
251                 Exception: if device not connected or data acquisition init failed.
252         """
253         self.__check_dev_status()
254         self.__with_q = with_q
255         if self.__status == iRecorder.Dev.SIGNAL:
256             return
257         if self.__status == iRecorder.Dev.IMPEDANCE:
258             self.stop_acquisition()
259         self.__status = iRecorder.Dev.SIGNAL_START
260         while self.__status not in [iRecorder.Dev.SIGNAL, iRecorder.Dev.TERMINATE]:
261             time.sleep(0.01)
262             self.__check_dev_status()
263

```

```

264     def get_data(
265         self, timeout: Optional[float] = 0.02
266     ) -> Optional[list[Optional[list]]]:
267         """
268             Acquire all available data, make sure this function is called in a loop when `with_q` is set to `True` in `start_acquisition_data()`
269
270         Args:
271             timeout: Non-negative value, blocks at most `timeout` seconds and return, if set to `None`, blocks until new data is available.
272
273         Returns:
274             A list of frames, each frame is a list contains all wanted eeg channels and trigger box channel,
275                 eeg channels can be updated by `update_channels()` .
276
277         Data Unit:
278             - eeg: micro volts ( $\mu$ V)
279             - triggerbox: int, from `0` to `255` .
280
281         Raises:
282             Exception: if device not connected or in data acquisition mode.
283
284             self.__check_dev_status()
285             if not self.__with_q:
286                 return
287             if self.__status != iRecorder.Dev.SIGNAL:
288                 raise Exception("Data acquisition not started, please start first.")
289             try:
290                 data: list = self.__save_data.get(timeout=timeout)
291             except queue.Empty:
292                 return []
293             while not self.__save_data.empty():
294                 data.extend(self.__save_data.get())
295             return data
296
297         def stop_acquisition(self) -> None:
298             """
299             Stop data or impedance acquisition, block until data acquisition stopped or failed.
300
301             Raises:
302                 Exception: if device not connected or acquisition stop failed.
303
304                 self.__check_dev_status()
305                 if self.__status == iRecorder.Dev.IDLE:
306                     return
307                 self.__status = iRecorder.Dev.IDLE_START
308                 while self.__status not in [iRecorder.Dev.IDLE, iRecorder.Dev.TERMINATE]:
309                     time.sleep(0.01)
310                 self.__check_dev_status()
311
312         def start_acquisition_impedance(self) -> None:
313             """
314             Send impedance acquisition command to device, block until data acquisition started or failed.
315
316             Raises:
317                 Exception: if device not connected or impedance acquisition init failed.
318
319                 self.__check_dev_status()
320                 if self.__status == iRecorder.Dev.IMPEDANCE:
321                     return
322                 if self.__status == iRecorder.Dev.SIGNAL:
323
324
325
326

```

```

327         self.stop_acquisition()
328     self.__status = iRecorder.Dev.IMPEDANCE_START
329     while self.__status not in [iRecorder.Dev.IMPEDANCE, iRecorder.Dev.TERMINATE]:
330         time.sleep(0.01)
331         self.__check_dev_status()
332
333     def get_impedance(self) -> Optional[list]:
334         """
335             Acquire channel impedances, return immediately, impedance update interval is about
336             2000ms.
337
338             Returns:
339                 A list of channel impedance ranging from `0` to `np.inf` if available, otherwise
340                 `None`.
341
342             Data Unit:
343                 - impedance: ohm ( $\Omega$ )
344         """
345         self.__check_dev_status()
346         return self.__parser.impedance
347
348     def close_dev(self) -> None:
349         """
350             Close device connection and release resources, resources are automatically released
351             on device error.
352         """
353         if self.__status != iRecorder.Dev.TERMINATE:
354             # ensure socket is closed correctly
355             self.__status = iRecorder.Dev.TERMINATE_START
356             while self.__status != iRecorder.Dev.TERMINATE:
357                 time.sleep(0.01)
358             if self.is_alive():
359                 self.join()
360
361     def get_packet_drop_times(self) -> int:
362         """
363             Retrieve packet drop times.
364             This value accumulates during data transmission and will be reset to `0` after device
365             status change.
366
367             Returns:
368                 accumulated packet drop times.
369         """
370         return self.__parser._drop_count
371
372     def get_battery_value(self) -> int:
373         """
374             Query battery level.
375
376             Returns:
377                 battery level in percentage, range from `0` to `100`.
378         """
379         return self.__parser.batt_val
380
381     def open_lsl_stream(self):
382         """
383             Open LSL stream, can be invoked after `start_acquisition_data()`,
384             each frame is the same as described in `get_data()` .
385
386             Raises:
387                 Exception: if data acquisition not started or LSL stream already opened.
388                 LSLError: if LSL stream creation failed.
389                 ImportError: if `pylsl` not installed or liblsl not installed on unix like

```

```

390     system.
391     """
392     if self.__status != iRecorder.Dev.SIGNAL:
393         raise Exception("Data acquisition not started, please start first.")
394     if hasattr(self, "_lsl_stream"):
395         raise Exception("LSL stream already opened.")
396     from ..utils.lslWrapper import lslSender
397
398     self._lsl_stream = lslSender(
399         self.__dev_args["ch_info"],
400         f"iRe{self.__dev_args['type']}_{self.__dev_args['name'][-2:]}",
401         "EEG",
402         self.__dev_args["fs"],
403         with_trigger=True,
404     )
405     self.__lsl_flag = True
406
407     def close_lsl_stream(self):
408         """
409             Close LSL stream manually, invoked automatically after `stop_acquisition()` or
410             `close_dev()`.
411         """
412         self.__lsl_flag = False
413         if hasattr(self, "_lsl_stream"):
414             del self._lsl_stream
415
416     def create_bdf_file(self, filename: str):
417         """
418             Create a BDF file and save data to it, invoke it after `start_acquisition_data()`.

419             Args:
420                 filename: file name to save data, accept absolute or relative path.
421
422             Raises:
423                 Exception: if data acquisition not started or `save_bdf_file` is invoked and BDF
424                 file already created.
425                 OSError: if BDF file creation failed, this may be caused by invalid file path or
426                 permission issue.
427                 ImportError: if `pyedflib` is not installed.
428
429             if self.__status != iRecorder.Dev.SIGNAL:
430                 raise Exception("Data acquisition not started")
431             if self._bdf_file is not None:
432                 raise Exception("BDF file already created.")
433             from ..utils.bdfWrapper import bdfSaverIRecorder
434
435             if filename[-4:].lower() != ".bdf":
436                 filename += ".bdf"
437             self._bdf_file = bdfSaverIRecorder(
438                 filename,
439                 self.__dev_args["ch_info"],
440                 self.__dev_args["fs"],
441                 f"iRecorder_{self.__dev_args['type']}_{self.__dev_args['name']}",
442             )
443             self.__bdf_flag = True
444
445     def close_bdf_file(self):
446         """
447             Close and save BDF file manually, invoked automatically after `stop_acquisition()` or
448             `close_dev()`.
449         """
450             self.__bdf_flag = False
451             if self._bdf_file is not None:

```

```
453         self._bdf_file.close_bdf()
454         self._bdf_file = None
455
456     def send_bdf_marker(self, marker: str):
457         """
458             Send marker to BDF file, can be invoked after `create_bdf_file()`, otherwise it will
459             be ignored.
460
461             Args:
462                 marker: marker string to write.
463
464                 if self._bdf_file is not None:
465                     self._bdf_file.write_annotation(marker)
466
467             # def set_callback_handler(self, handler: Callable[[Optional[str]], None]):
468             #     """
469             #     Set callback handler function, invoked automatically when device thread ended if
470             #     set.
471
472             #     Args:
473             #         handler: a callable function that takes a string of error message or `None` as
474             #     input.
475             #         """
476             #         self.handler = handler
477
478         def __check_dev_status(self):
479             if self.__error_message is None:
480                 return
481             if self.is_alive():
482                 self.close_dev()
483             raise Exception(self.__error_message)
484
485     def run(self):
486         while self.__status not in [iRecorder.Dev.TERMINATE_START]:
487             if self.__status == iRecorder.Dev.SIGNAL_START:
488                 self.__recv_data(imp_mode=False)
489             elif self.__status == iRecorder.Dev.IMPEDANCE_START:
490                 self.__recv_data(imp_mode=True)
491             elif self.__status in [iRecorder.Dev.IDLE_START]:
492                 self.__idle_state()
493             else:
494                 self.__error_message = f"Unknown status: {self.__status}"
495                 break
496             try:
497                 self.dev.close_socket()
498             except Exception:
499                 pass
500             finally:
501                 self.__finish_search()
502                 self.__status = iRecorder.Dev.TERMINATE
503                 # if self.handler is not None:
504                 #     self.handler(self.__error_message)
505
506     def __recv_data(self, imp_mode=True):
507         self.__parser.imp_flag = imp_mode
508         retry = 0
509         try:
510             if imp_mode:
511                 self.dev.start_impe()
512                 self.__status = iRecorder.Dev.IMPEDANCE
513                 print("IMPEDANCE START")
514             else:
515                 self.dev.start_data()
```

```

516         self.__status = iRecorder.Dev.SIGNAL
517         print("SIGNAL START")
518     except Exception:
519         self.__error_message = "Data/Impedance mode initialization failed."
520         self.__status = iRecorder.Dev.TERMINATE_START
521     # recv data
522     while self.__status in [iRecorder.Dev.SIGNAL, iRecorder.Dev.IMPEDANCE]:
523         try:
524             data = self.dev.recv_socket()
525             if not data:
526                 raise Exception("Remote end closed.")
527             ret = self.__parser.parse_data(data)
528             if ret:
529                 if self.__with_q:
530                     self.__save_data.put(ret)
531                 if self.__bdf_flag:
532                     self._bdf_file.write_chunk(ret)
533                 if self.__lsl_flag:
534                     self._lsl_stream.push_chunk(ret)
535         except Exception:
536             if (self.__dev_args["type"] == "W32") and (retry < 1):
537                 try:
538                     print("Wi-Fi reconnecting...")
539                     self.dev.close_socket()
540                     self.dev = self.__dev_sock(self.__dev_args, retry_timeout=3)
541                     retry += 1
542                     continue
543                 except Exception:
544                     print("Wi-Fi reconnection failed")
545                     self.__error_message = "Data transmission timeout."
546                     self.__status = iRecorder.Dev.TERMINATE_START
547             # postprocess
548             self.close_bdf_file()
549             self.close_lsl_stream()
550             self.__parser.clear_buffer()
551             while not self.__save_data.empty():
552                 self.__save_data.get()
553             # stop recv data
554             if self.__status != iRecorder.Dev.TERMINATE_START:
555                 try: # stop data acquisition when thread ended
556                     self.dev.stop_recv()
557                 except Exception:
558                     if self.__status == iRecorder.Dev.IDLE_START:
559                         self.__error_message = "Device connection lost."
560                     self.__status = iRecorder.Dev.TERMINATE_START

def __idle_state(self):
    timestamp = time.time()
    self.__status = iRecorder.Dev.IDLE
    while self.__status in [iRecorder.Dev.IDLE]:
        if (time.time() - timestamp) < 5:
            time.sleep(0.2) # to reduce cpu usage
            continue
        try: # heartbeat to keep socket alive and update battery level
            self.__parser.batt_val = self.dev.send_heartbeat()
            timestamp = time.time()
            # print("Ah, ah, ah, ah\nStayin' alive, stayin' alive")
        except Exception:
            self.__error_message = "Device connection lost."
            self.__status = iRecorder.Dev.TERMINATE_START

def __get_chs(self) -> int:
    return int("".join([i for i in self.__dev_args["type"] if i.isdigit()]))

```

```
def __finish_search(self):
    if self.__interface.is_alive():
        self.__interface.stop()
        self.__interface.join()
    while not self.__info_q.empty():
        self.__info_q.get()
```

2.1.1 __init__(dev_type)

Parameters:

Name	Type	Description	Default
dev_type	str	iRecorder device type. available options: Literal["W8", "USB8", "W16", "USB16", "W32", "USB32"]	required

Raises: Exception: if device type not supported. Exception: if adapter not available.

Source code in `src\iRecorder\device.py`

```
24     def __init__(  
25         self,  
26         dev_type: str  
27     ):  
28         """  
29         Args:  
30             dev_type: iRecorder device type. available options: Literal["W8", "USB8", "W16",  
31             "USB16", "W32", "USB32"]  
32         Raises:  
33             Exception: if device type not supported.  
34             Exception: if adapter not available.  
35         """  
36         if dev_type not in {"W8", "USB8", "W16", "USB16", "W32", "USB32"}:  
37             raise ValueError("Unsupported device type.")  
38         super().__init__(daemon=True, name=f"iRecorder {dev_type}")  
39         self.handler = None  
40         self.__info_q = Queue(128)  
41         self.__with_q = True  
42         self.__error_message = "Device not connected, please connect first."  
43         self.__save_data = Queue()  
44         self.__status = iRecorder.Dev.TERMINATE  
45         self.__lsl_flag = False  
46         self.__bdf_flag = False  
47         self.__dev_args = {"type": dev_type}  
48         self.__dev_args.update({"channel": self.__get_chs()})  
49  
50         self.__parser = Parser(self.__dev_args["channel"])  
51         self.__interface = get_interface(dev_type)(self.__info_q)  
52         self.__dev_sock = get_sock(dev_type)  
53         self.__dev_args.update({"AdapterInfo": self.__interface.interface})  
54  
55         self._bdf_file = None  
56         self.dev = None  
57  
58         self.set_frequency()  
      self.update_channels()
```

2.1.2 `close_bdf_file()`

Close and save BDF file manually, invoked automatically after `stop_acquisition()` or `close_dev()`

Source code in `src\iRecorder\device.py`

```
420     def close_bdf_file(self):  
421         """  
422             Close and save BDF file manually, invoked automatically after `stop_acquisition()` or  
423             `close_dev()`  
424         """  
425             self.__bdf_flag = False  
426             if self._bdf_file is not None:  
427                 self._bdf_file.close_bdf()  
                 self._bdf_file = None
```

2.1.3 close_dev()

Close device connection and release resources, resources are automatically released on device error.

Source code in `src\iRecorder\device.py`

```

328     def close_dev(self) -> None:
329         """
330             Close device connection and release resources, resources are automatically released on
331             device error.
332         """
333         if self.__status != iRecorder.Dev.TERMINATE:
334             # ensure socket is closed correctly
335             self.__status = iRecorder.Dev.TERMINATE_START
336             while self.__status != iRecorder.Dev.TERMINATE:
337                 time.sleep(0.01)
338         if self.is_alive():
339             self.join()

```

2.1.4 close_lsl_stream()

Close LSL stream manually, invoked automatically after `stop_acquisition()` or `close_dev()`

Source code in `src\iRecorder\device.py`

```

384     def close_lsl_stream(self):
385         """
386             Close LSL stream manually, invoked automatically after `stop_acquisition()` or
387             `close_dev()`.
388         """
389         self.__lsl_flag = False
390         if hasattr(self, "_lsl_stream"):
391             del self._lsl_stream

```

2.1.5 connect_device(addr)

Connect to device by address, block until connection is established or failed.

Parameters:

Name	Type	Description	Default
addr	str	device address.	<i>required</i>

Raises:

Type	Description
Exception	if device already connected or connection establishment failed.

Source code in <code>src\iRecorder\device.py</code>	
	<pre> 180 def connect_device(self, addr: str) -> None: 181 """ 182 Connect to device by address, block until connection is established or failed. 183 184 Args: 185 addr: device address. 186 187 Raises: 188 Exception: if device already connected or connection establishment failed. 189 """ 190 if self.is_alive(): 191 raise Exception("iRecorder already connected") 192 try: 193 ret = self.__interface.connect(addr) 194 self.__dev_args.update({"name": addr, "sock": ret}) 195 self.__dev_args.update({"_length": self.__parser.packet_len}) 196 self.dev = self.__dev_sock(self.__dev_args) 197 self.__parser.batt_val = self.dev.send_heartbeat() 198 self.__error_message = None 199 self.__status = iRecorder.Dev.IDLE_START 200 self.start() 201 except Exception as e: 202 self.__error_message = "Device connection failed." 203 self.__finish_search() 204 raise e </pre>

2.1.6 `create_bdf_file(filename)`

Create a BDF file and save data to it, invoke it after `start_acquisition_data()`.

Parameters:

Name	Type	Description	Default
filename	str	file name to save data, accept absolute or relative path.	required

Raises:

Type	Description
Exception	if data acquisition not started or <code>save_bdf_file</code> is invoked and BDF file already created.

Type	Description
OSError	if BDF file creation failed, this may be caused by invalid file path or permission issue.
importError	if <code>pyedflib</code> is not installed.

Source code in `src\ConEXG\iRecorder\device.py`

```

392     def create_bdf_file(self, filename: str):
393         """
394             Create a BDF file and save data to it, invoke it after `start_acquisition_data()`.
395
396             Args:
397                 filename: file name to save data, accept absolute or relative path.
398
399             Raises:
400                 Exception: if data acquisition not started or `save_bdf_file` is invoked and BDF file
401 already created.
402                 OSError: if BDF file creation failed, this may be caused by invalid file path or
403 permission issue.
404                 ImportError: if `pyedflib` is not installed.
405
406             if self.__status != iRecorder.Dev.SIGNAL:
407                 raise Exception("Data acquisition not started")
408             if self._bdf_file is not None:
409                 raise Exception("BDF file already created.")
410             from ..utils.bdfWrapper import bdfSaverIRecorder
411
412             if filename[-4:].lower() != ".bdf":
413                 filename += ".bdf"
414             self._bdf_file = bdfSaverIRecorder(
415                 filename,
416                 self.__dev_args["ch_info"],
417                 self.__dev_args["fs"],
418                 f"iRecorder_{self.__dev_args['type']}_{self.__dev_args['name']}",
419             )
420             self.__bdf_flag = True

```

2.1.7 `find_devs(duration=None)`

Search for available devices, can only be called once per instance.

Parameters:

Name	Type	Description	Default
duration	Optional[int]	Search interval in seconds, blocks for about <code>duration</code> seconds and return found devices, if set to <code>None</code> , return	None

Name	Type	Description	Default
	None	immediately, devices can later be acquired by calling <code>get_devs()</code> in a loop.	

Returns:

Type	Description
Optional[list]	Available devices.

Raises:

Type	Description
Exception	If search thread already running or iRecorder already connected.

Source code in `src\iRecorder\device.py`

```

60 def find_devs(self, duration: Optional[int] = None) -> Optional[list]:
61     """
62         Search for available devices, can only be called once per instance.
63
64     Args:
65         duration: Search interval in seconds, blocks for about `duration` seconds and return
66         found devices,
67             if set to `None`, return `None` immediately, devices can later be acquired by
68         calling `get_devs()` in a loop.
69
70     Returns:
71         Available devices.
72
73     Raises:
74         Exception: If search thread already running or iRecorder already connected.
75     """
76     if self.is_alive():
77         raise Exception("iRecorder already connected.")
78     if self.__interface.is_alive():
79         raise Exception("Search thread already running.")
80     self.__interface.start()
81     if duration is None:
82         return
83     start = time.time()
84     while time.time() - start < duration:
85         time.sleep(0.5)
86     self.__finish_search()
87     return self.get_devs()

```

2.1.8 `get_available_frequency(dev_type)` staticmethod

Get available sample frequencies of different device types.

Returns:

Type	Description
list	Available sample frequencies in Hz.

Source code in `src\ConEXG\iRecorder\device.py`

```

141     @staticmethod
142     def get_available_frequency(dev_type: str) -> list:
143         """Get available sample frequencies of different device types.
144
145         Returns:
146             Available sample frequencies in Hz.
147         """
148         if "USB" in dev_type:
149             return [500, 1000, 2000]
150         return [500]
```

2.1.9 `get_battery_value()`

Query battery level.

Returns:

Type	Description
int	battery level in percentage, range from 0 to 100.

Source code in `src\ConEXG\iRecorder\device.py`

```

350     def get_battery_value(self) -> int:
351         """
352             Query battery level.
353
354         Returns:
355             battery level in percentage, range from `0` to `100`.
356         """
357         return self.__parser.batt_val
```

2.1.10 `get_data(timeout=0.02)`

Acquire all available data, make sure this function is called in a loop when `with_q` is set to `True` in `start_acquisition_data()`

Parameters:

Name	Type	Description	Default
<code>timeout</code>	<code>Optional[float]</code>	Non-negative value, blocks at most <code>timeout</code> seconds and return, if set to <code>None</code> , blocks until new data is available.	<code>0.02</code>

Returns:

Type	Description
<code>Optional[list[Optional[list]]]</code>	A list of frames, each frame is a list contains all wanted eeg channels and trigger box channel, eeg channels can be updated by <code>update_channels()</code> .

Data Unit ▼

- eeg: micro volts (μ V)
- triggerbox: int, from 0 to 255

Raises:

Type	Description
<code>Exception</code>	if device not connected or in data acquisition mode.

Source code in `src\iRecorder\device.py`

```

250     def get_data(
251         self, timeout: Optional[float] = 0.02
252     ) -> Optional[list[Optional[list]]]:
253         """
254             Acquire all available data, make sure this function is called in a loop when `with_q` is
255             set to `True` in `start_acquisition_data()`.
256
257             Args:
258                 timeout: Non-negative value, blocks at most `timeout` seconds and return, if set to
259                 `None`, blocks until new data is available.
260
261             Returns:
262                 A list of frames, each frame is a list contains all wanted eeg channels and trigger
263                 box channel,
264                     eeg channels can be updated by `update_channels()`.
265
266             Data Unit:
267                 - eeg: micro volts ( $\mu$ V)
268                 - triggerbox: int, from `0` to `255`
269
270             Raises:
271                 Exception: if device not connected or in data acquisition mode.
272             """
273             self.__check_dev_status()
274             if not self.__with_q:
275                 return
276             if self.__status != iRecorder.Dev.SIGNAL:
277                 raise Exception("Data acquisition not started, please start first.")
278             try:
279                 data: list = self.__save_data.get(timeout=timeout)
280             except queue.Empty:
281                 return []
282             while not self.__save_data.empty():
283                 data.extend(self.__save_data.get())
284             return data

```

2.1.11 `get_dev_info()`

Get current device information, including device name, hardware channel number, acquired channels, sample frequency, etc.

Returns:

Type	Description
dict	A dictionary containing device information, which includes: <code>type</code> : hardware type; <code>channel</code> : hardware channel number; <code>AdapterInfo</code> : adapter used for connection; <code>fs</code> : sample frequency in Hz; <code>ch_info</code> : channel dictionary, including channel index and name, can be altered by <code>update_channels()</code> .

Source code in `src\iRecorder\device.py`

```

127 def get_dev_info(self) -> dict:
128     """
129         Get current device information, including device name, hardware channel number, acquired
130         channels, sample frequency, etc.
131
132     Returns:
133         A dictionary containing device information, which includes:
134             `type`: hardware type;
135             `channel`: hardware channel number;
136             `AdapterInfo`: adapter used for connection;
137             `fs`: sample frequency in Hz;
138             `ch_info`: channel dictionary, including channel index and name, can be altered
139         by `update_channels()``.
140
141     """
142
143     return deepcopy(self.__dev_args)

```

2.1.12 `get_dev_status()`

Get current device status.

Returns:

Type	Description
str	"SIGNAL": data acquisition mode
str	"IMPEDANCE": impedance acquisition mode
str	"IDLE": idle mode
str	"TERMINATE": device not connected or connection closed.

Source code in `src\iRecorder\device.py`

```

115 def get_dev_status(self) -> str:
116     """
117         Get current device status.
118
119     Returns:
120         "SIGNAL": data acquisition mode
121         "IMPEDANCE": impedance acquisition mode
122         "IDLE": idle mode
123         "TERMINATE": device not connected or connection closed.
124
125     """
126
127     return self.__status.name

```

2.1.13 `get_devs(verbose=False)`

Get available devices. This can be called after `find_devs(duration = None)` in a loop, each call will *only* return newly found devices.

Parameters:

Name	Type	Description	Default
verbose	bool	if True, return all available devices information, otherwise only return names for connection, if you don't know what this parameter does, just leave it at its default value.	False

Returns:

Type	Description
list	Newly found devices.

Raises:

Type	Description
Exception	adapter not found or not enabled etc.

Source code in `src\ eConEXG\ iRecorder\ device.py`

```

87 def get_devs(self, verbose: bool = False) -> list:
88     """
89     Get available devices. This can be called after `find_devs(duration = None)` in a loop,
90     each call will *only* return newly found devices.
91
92     Args:
93         verbose: if True, return all available devices information, otherwise only return
94         names for connection,
95             if you don't know what this parameter does, just leave it at its default value.
96
97     Returns:
98         Newly found devices.
99
100    Raises:
101        Exception: adapter not found or not enabled etc.
102    """
103    ret = []
104    time.sleep(0.1)
105    while not self.__info_q.empty():
106        info = self.__info_q.get()
107        if isinstance(info, list):
108            ret.append(info if verbose else info[-1])
109        elif isinstance(info, bool):
110            if verbose:
111                ret.append(info)
112            elif isinstance(info, str):
113                raise Exception(info)
114    return ret

```

2.1.14 `get_impedance()`

Acquire channel impedances, return immediately, impedance update interval is about 2000ms.

Returns:

Type	Description
<code>Optional[list]</code>	A list of channel impedance ranging from <code>0</code> to <code>np.inf</code> if available, otherwise <code>None</code> .

Data Unit

- `impedance: ohm (Ω)`

Source code in `src\iRecorder\device.py`

```

315 def get_impedance(self) -> Optional[list]:
316     """
317         Acquire channel impedances, return immediately, impedance update interval is about
318         2000ms.
319
320     Returns:
321         A list of channel impedance ranging from `0` to `np.inf` if available, otherwise
322         `None`.
323
324     Data Unit:
325         - impedance: ohm ( $\Omega$ )
326         """
327     self.__check_dev_status()
328     return self.__parser.impedance

```

`2.1.15 get_packet_drop_times()`

Retrieve packet drop times. This value accumulates during data transmission and will be reset to `0` after device status change.

Returns:

Type	Description
<code>int</code>	accumulated packet drop times.

Source code in `src\iRecorder\device.py`

```

340 def get_packet_drop_times(self) -> int:
341     """
342         Retrieve packet drop times.
343         This value accumulates during data transmission and will be reset to `0` after device
344         status change.
345
346     Returns:
347         accumulated packet drop times.
348         """
349     return self.__parser._drop_count

```

`2.1.16 open_lsl_stream()`

Open LSL stream, can be invoked after `start_acquisition_data()`, each frame is the same as described in `get_data()`.

Raises:

Type	Description
Exception	if data acquisition not started or LSL stream already opened.
LSLException	if LSL stream creation failed.
importError	if `pylsl` not installed or `liblsl` not installed on unix like system.

Source code in `src\iRecorder\device.py`

```

359 def open_lsl_stream(self):
360     """
361     Open LSL stream, can be invoked after `start_acquisition_data()`,
362     each frame is the same as described in `get_data()`.
363
364     Raises:
365         Exception: if data acquisition not started or LSL stream already opened.
366         LSLException: if LSL stream creation failed.
367         ImportError: if `pylsl` not installed or `liblsl` not installed on unix like system.
368     """
369     if self.__status != iRecorder.Dev.SIGNAL:
370         raise Exception("Data acquisition not started, please start first.")
371     if hasattr(self, "_lsl_stream"):
372         raise Exception("LSL stream already opened.")
373     from ..utils.lslWrapper import lslSender
374
375     self._lsl_stream = lslSender(
376         self.__dev_args["ch_info"],
377         f"iRe{self.__dev_args['type']}_{self.__dev_args['name'][-2:]}",
378         "EEG",
379         self.__dev_args["fs"],
380         with_trigger=True,
381     )
382     self.__lsl_flag = True

```

2.1.17 `send_bdf_marker(marker)`

Send marker to BDF file, can be invoked after `create_bdf_file()`, otherwise it will be ignored.

Parameters:

Name	Type	Description	Default
marker	str	marker string to write.	<i>required</i>

Source code in `src\iRecorder\device.py`

```

429 def send_bdf_marker(self, marker: str):
430     """
431         Send marker to BDF file, can be invoked after `create_bdf_file()`, otherwise it will be
432         ignored.
433
434     Args:
435         marker: marker string to write.
436     """
437     if self._bdf_file is not None:
438         self._bdf_file.write_annotation(marker)

```

2.1.18 `set_frequency(fs=None)`

Update device sample frequency, this method should be invoked before `connect_device`.

Parameters:

Name	Type	Description	Default
<code>fs</code>	<code>Optional[int]</code>	sample frequency in Hz, if <code>fs</code> is set to <code>None</code> or not in <code>get_available_frequency()</code> , it will fall back to the lowest available frequency.	<code>None</code>

Raises:

Type	Description
<code>Exception</code>	Device is already connected.

 New in

- now you can set the sample frequency after device connection.

Source code in `src\iRecorder\device.py`

```

152     def set_frequency(self, fs: Optional[int] = None):
153         """Update device sample frequency, this method should be invoked before `connect_device`.
154
155         Args:
156             fs: sample frequency in Hz, if `fs` is set to `None` or not in
157             `get_available_frequency()`,
158                 it will fall back to the lowest available frequency.
159
160         Raises:
161             Exception: Device is already connected.
162
163         New in:
164             - now you can set the sample frequency after device connection.
165         """
166         if self.__status not in [iRecorder.Dev.IDLE, iRecorder.Dev.TERMINATE]:
167             warn = "Device acquisition in progress, please `stop_acquisition()` first."
168             raise Exception(warn)
169         available = self.get_available_frequency(self.__dev_args["type"])
170         default = available[0]
171         if fs is None:
172             fs = default
173         if fs not in available:
174             print(f"Invalid sample frequency, fallback to {default}Hz")
175             fs = default
176         self.__dev_args.update({"fs": fs})
177         self.__parser._update_fs(fs)
178         if self.dev is not None:
179             self.dev.set_fs(fs)

```

2.1.19 `start_acquisition_data(with_q=True)`

Send data acquisition command to device, block until data acquisition started or failed.

Parameters:

Name	Type	Description	Default
<code>with_q</code>	<code>bool</code>	if True, signal data will be stored in a queue and should be acquired by calling <code>get_data()</code> in a loop in case data queue is full. if False, new data will not be directly available and can only be acquired through <code>open_lsl_stream</code> and <code>save_bdf_file</code> .	<code>True</code>

Raises:

Type	Description
<code>Exception</code>	if device not connected or data acquisition init failed.

Source code in [src\iRecorder\device.py](#)

```
228 def start_acquisition_data(self, with_q: bool = True):
229     """
230         Send data acquisition command to device, block until data acquisition started or failed.
231
232     Args:
233         with_q: if True, signal data will be stored in a queue and **should** be acquired by
234         calling `get_data()` in a loop in case data queue is full.
235             if False, new data will not be directly available and can only be acquired
236         through `open_lsl_stream` and `save_bdf_file`.
237
238     Raises:
239         Exception: if device not connected or data acquisition init failed.
240     """
241     self.__check_dev_status()
242     self.__with_q = with_q
243     if self.__status == iRecorder.Dev.SIGNAL:
244         return
245     if self.__status == iRecorder.Dev.IMPEDANCE:
246         self.stop_acquisition()
247     self.__status = iRecorder.Dev.SIGNAL_START
248     while self.__status not in [iRecorder.Dev.SIGNAL, iRecorder.Dev.TERMINATE]:
249         time.sleep(0.01)
250     self.__check_dev_status()
```

2.1.20 `start_acquisition_impedance()`

Send impedance acquisition command to device, block until data acquisition started or failed.

Raises:

Type	Description
Exception	if device not connected or impedance acquisition init failed.

” Source code in `src\iRecorder\device.py`

```

298 def start_acquisition_impedance(self) -> None:
299     """
300         Send impedance acquisition command to device, block until data acquisition started or
301         failed.
302
303     Raises:
304         Exception: if device not connected or impedance acquisition init failed.
305     """
306     self.__check_dev_status()
307     if self.__status == iRecorder.Dev.IMPEDANCE:
308         return
309     if self.__status == iRecorder.Dev.SIGNAL:
310         self.stop_acquisition()
311     self.__status = iRecorder.Dev.IMPEDANCE_START
312     while self.__status not in [iRecorder.Dev.IMPEDANCE, iRecorder.Dev.TERMINATE]:
313         time.sleep(0.01)
314     self.__check_dev_status()

```

2.1.21 `stop_acquisition()`

Stop data or impedance acquisition, block until data acquisition stopped or failed.

Raises:

Type	Description
Exception	if device not connected or acquisition stop failed.

” Source code in `src\iRecorder\device.py`

```

283 def stop_acquisition(self) -> None:
284     """
285         Stop data or impedance acquisition, block until data acquisition stopped or failed.
286
287     Raises:
288         Exception: if device not connected or acquisition stop failed.
289     """
290     self.__check_dev_status()
291     if self.__status == iRecorder.Dev.IDLE:
292         return
293     self.__status = iRecorder.Dev.IDLE_START
294     while self.__status not in [iRecorder.Dev.IDLE, iRecorder.Dev.TERMINATE]:
295         time.sleep(0.01)
296     self.__check_dev_status()

```

2.1.22 `update_channels(channels=None)`

Update channels to acquire, invoke this method when device is not acquiring data or impedance.

Parameters:

Name	Type	Description	Default
channels	Optional[dict]	channel number and name mapping, e.g. {0: "FPz", 1: "Oz", 2: "CPz"}, if None is given, reset to all available channels with default names.	None

Raises:

Type	Description
Exception	if data/impedance acquisition in progress.

Source code in [src\iRecorder\device.py](#)

```

206     def update_channels(self, channels: Optional[dict] = None):
207         """
208             Update channels to acquire, invoke this method when device is not acquiring data or
209             impedance.
210
211             Args:
212                 channels: channel number and name mapping, e.g. `{'0': 'FPz', 1: 'Oz', 2: 'CPz'}`,
213                     if `None` is given, reset to all available channels with default names.
214
215             Raises:
216                 Exception: if data/impedance acquisition in progress.
217             """
218             if self.__status not in [iRecorder.Dev.IDLE, iRecorder.Dev.TERMINATE]:
219                 warn = "Device acquisition in progress, please stop_acquisition() first."
220                 raise Exception(warn)
221             if channels is None:
222                 from .default_config import getChannels
223
224                 channels = getChannels(self.__dev_args["channel"])
225                 self.__dev_args.update({"ch_info": channels})
226                 ch_idx = [i for i in channels.keys()]
227                 self.__parser._update_chs(ch_idx)

```

3 iFocus

3.1 iFocus

Bases: Thread

Source code in `src\ConEXG\iFocus\data_reader.py`

```

13  class iFocus(Thread):
14      class Dev(Enum):
15          SIGNAL = 10
16          SIGNAL_START = 11
17          IDLE = 30
18          IDLE_START = 31
19          TERMINATE = 40
20          TERMINATE_START = 41
21
22      dev_args = {
23          "type": "iFocus",
24          "fs_eeg": 250,
25          "fs_imu": 50,
26          "channel_eeg": {0: "CH0"},  

27          "channel_imu": {0: "X", 1: "Y", 2: "Z"},  

28          "AdapterInfo": "Serial Port",
29      }
30
31  def __init__(self, port: Optional[str] = None) -> None:
32      """
33          Args:
34              port: if not given, connect to the first available device.
35      """
36      super().__init__(daemon=True)
37      self.__status = iFocus.Dev.TERMINATE
38      if port is None:
39          port = iFocus.find_devs()[0]
40      self.__save_data = Queue()
41      self.__parser = Parser()
42      self.dev_args = deepcopy(iFocus.dev_args)
43      self.dev = sock(port)
44      self.set_frequency()
45      self.__with_q = True
46      self.__socket_flag = "Device not connected, please connect first."
47      self.__bdf_flag = False
48      try:
49          self.dev.connect_socket()
50      except Exception as e:
51          try:
52              self.dev.close_socket()
53          finally:
54              raise e
55      self.__status = iFocus.Dev.IDLE_START
56      self.__socket_flag = None
57      self._lsl_eeg = None
58      self._lsl_imu = None
59      self.__lsl_imu_flag = False
60      self.__lsl_eeg_flag = False
61      self._bdf_file = None
62      self.__enable_imu = False
63      self.dev_args["name"] = port
64      self.start()
65
66  def set_frequency(self, fs_eeg: int = None):
67      """
68          Change the sampling frequency of iFocus.
69
70          Args:
71              fs_eeg: sampling frequency of eeg data, should be 250 or 500,
72                      fs_imu will be automatically set to 1/5 of fs_eeg.
73
74          Raises:

```

```

75             ValueError: if fs_eeg is not 250 or 500.
76             NotImplementedError: device firmware too old, not supporting 500Hz.
77         """
78         if self.__status == iFocus.Dev.SIGNAL:
79             raise Exception("Data acquisition already started, please stop first.")
80         if fs_eeg is None:
81             fs_eeg = self.dev_args["fs_eeg"]
82         if fs_eeg not in [250, 500]:
83             raise ValueError("fs_eeg should be 250 or 500")
84         self.dev_args["fs_eeg"] = fs_eeg
85         fs_imu = fs_eeg // 5
86         self.dev_args["fs_imu"] = fs_imu
87         if hasattr(self, "dev"):
88             self.dev.set_frequency(fs_eeg)
89
90     def get_dev_info(self) -> dict:
91         """
92             Get current device information, including device name, hardware channel number,
93             acquired channels, sample frequency, etc.
94
95             Returns:
96                 A dictionary containing device information, which includes:
97                     `type`: hardware type;
98                     `channel_eeg`: channel dictionary, including EEG channel index and name;
99                     `channel_imu`: channel dictionary, including IMU channel index and name;
100                    `AdapterInfo`: adapter used for connection;
101                    `fs_eeg`: sample frequency of EEG in Hz;
102                    `fs_imu`: sample frequency of IMU in Hz;
103
104            return deepcopy(self.dev_args)
105
106    @staticmethod
107    def find_devs() -> list:
108        """
109            Find available iFocus devices.
110
111            Returns:
112                available device ports.
113
114            Raises:
115                Exception: if no iFocus device found.
116
117            return sock._find_devs()
118
119    def get_data(
120        self, timeout: Optional[float] = 0.02
121    ) -> Optional[list[Optional[list]]]:
122        """
123            Acquire all available data, make sure this function is called in a loop when `with_q`
124            is set to `True` in `start_acquisition_data()`
125
126            Args:
127                timeout: Non-negative value, blocks at most 'timeout' seconds and return, if set
128                to `None`, blocks until new data available.
129
130            Returns:
131                A list of frames, each frame is made up of 5 eeg data and 1 imu data in a shape
132                as below:
133                [[`eeg_0`], [`eeg_1`], [`eeg_2`], [`eeg_3`], [`eeg_4`], [`imu_x`, `imu_y`,
134                `imu_z`]],
135                in which number `0~4` after `_` indicates the time order of channel data.
136
137            Raises:

```

```

138             Exception: if device not connected, connection failed, data transmission
139             timeout/init failed, or unknown error.
140
141         Data Unit:
142             - eeg: µV
143             - imu: degree(°)
144         """
145
146     self.__check_dev_status()
147     if not self.__with_q:
148         return
149     try:
150         data: list = self.__save_data.get(timeout=timeout)
151     except queue.Empty:
152         return []
153     while not self.__save_data.empty():
154         data.extend(self.__save_data.get())
155     return data
156
157     def start_acquisition_data(self, with_q: bool = True) -> None:
158         """
159             Send data acquisition command to device, block until data acquisition started or
160             failed.
161
162             Args:
163                 with_q: if True, signal data will be stored in a queue and **should** be acquired
164                 by calling `get_data()` in a loop in case data queue is full.
165                 if False, new data will not be directly available and can only be acquired
166                 through lsl stream.
167
168             """
169         self.__check_dev_status()
170         self.__with_q = with_q
171         if self.__status == iFocus.Dev.SIGNAL:
172             return
173         self.__status = iFocus.Dev.SIGNAL_START
174         while self.__status not in [iFocus.Dev.SIGNAL, iFocus.Dev.TERMINATE]:
175             time.sleep(0.01)
176         self.__check_dev_status()
177
178     def stop_acquisition(self) -> None:
179         """
180             Stop data or impedance acquisition, block until data acquisition stopped or failed.
181
182             """
183         self.__check_dev_status()
184         self.__status = iFocus.Dev.IDLE_START
185         while self.__status not in [iFocus.Dev.IDLE, iFocus.Dev.TERMINATE]:
186             time.sleep(0.01)
187         self.__check_dev_status()
188
189     def open_lsl_eeg(self):
190         """
191             Open LSL EEG stream, can be invoked after `start_acquisition_data()`.
192
193             Raises:
194                 Exception: if data acquisition not started or LSL stream already opened.
195                 LSLError: if LSL stream creation failed.
196                 ImportError: if `pylsl` is not installed or liblsl not installed for unix like
197                 system.
198
199                 if self.__status != iFocus.Dev.SIGNAL:
200                     raise Exception("Data acquisition not started, please start first.")
201                 if hasattr(self, "_lsl_eeg"):
202                     raise Exception("LSL stream already opened.")

```

```
201     from ..utils.lslWrapper import lslSender
202
203     self._lsl_eeg = lslSender(
204         self.dev_args["channel_eeg"],
205         f"{self.dev_args['type']}EEG{self.dev_args['name'][-2:]}",
206         "EEG",
207         self.dev_args["fs_eeg"],
208         with_trigger=False,
209     )
210     self.__lsl_eeg_flag = True
211
212     def close_lsl_eeg(self):
213         """
214             Close LSL EEG stream manually, invoked automatically after `stop_acquisition()` and
215             `close_dev()`
216         """
217         self.__lsl_eeg_flag = False
218         if hasattr(self, "_lsl_eeg"):
219             del self._lsl_eeg
220
221     def open_lsl_imu(self):
222         """
223             Open LSL IMU stream, can be invoked after `start_acquisition_data()` .
224
225             Raises:
226                 Exception: if data acquisition not started or LSL stream already opened.
227                 LSLError: if LSL stream creation failed.
228                 ImportError: if `pylsl` is not installed or liblsl not installed for unix like
229             system.
230         """
231         if self.__status != iFocus.Dev.SIGNAL:
232             raise Exception("Data acquisition not started, please start first.")
233         if hasattr(self, "_lsl_imu"):
234             raise Exception("LSL stream already opened.")
235         from ..utils.lslWrapper import lslSender
236
237         self._lsl_imu = lslSender(
238             self.dev_args["channel_imu"],
239             f"{self.dev_args['type']}IMU{self.dev_args['name'][-2:]}",
240             "IMU",
241             self.dev_args["fs_imu"],
242             unit="degree",
243             with_trigger=False,
244         )
245         self.__lsl_imu_flag = True
246
247     def close_lsl_imu(self):
248         """
249             Close LSL IMU stream manually, invoked automatically after `stop_acquisition()` and
250             `close_dev()`
251         """
252         self.__lsl_imu_flag = False
253         if hasattr(self, "_lsl_imu"):
254             del self._lsl_imu
255
256     def setIMUFlag(self, check):
257         self.__enable_imu = check
258
259     def create_bdf_file(self, filename: str):
260         """
261             Create a BDF file and save data to it, invoke it after `start_acquisition_data()` .
262
263             Args:

```

```

264         filename: file name to save data, accept absolute or relative path.
265
266     Raises:
267         Exception: if data acquisition not started or `save_bdf_file` is invoked and BDF
268         file already created.
269         OSError: if BDF file creation failed, this may be caused by invalid file path or
270         permission issue.
271         ImportError: if `pyedflib` is not installed.
272         """
273         if self.__status != iFocus.Dev.SIGNAL:
274             raise Exception("Data acquisition not started")
275         if self._bdf_file is not None:
276             raise Exception("BDF file already created.")
277         from ..utils.bdfWrapper import bdfSaverEEG, bdfSaverEEGIMU
278
279         if filename[-4:].lower() != ".bdf":
280             filename += ".bdf"
281         if self.__enable_imu:
282             self._bdf_file = bdfSaverEEGIMU(
283                 filename,
284                 self.dev_args["channel_eeg"],
285                 self.dev_args["fs_eeg"],
286                 self.dev_args["channel_imu"],
287                 self.dev_args["fs_imu"],
288                 self.dev_args["type"],
289             )
290         else:
291             self._bdf_file = bdfSaverEEG(
292                 filename,
293                 self.dev_args["channel_eeg"],
294                 self.dev_args["fs_eeg"],
295                 self.dev_args["type"],
296             )
297         self.__bdf_flag = True
298
299     def close_bdf_file(self):
300         """
301             Close and save BDF file manually, invoked automatically after `stop_acquisition()` or
302             `close_dev()`.
303         """
304         self.__bdf_flag = False
305         if self._bdf_file is not None:
306             self._bdf_file.close_bdf()
307             del self._bdf_file
308
309     def send_bdf_marker(self, marker: str):
310         """
311             Send marker to BDF file, can be invoked after `create_bdf_file()`, otherwise it will
312             be ignored.
313
314             Args:
315                 marker: marker string to write.
316             """
317             if hasattr(self, "_bdf_file"):
318                 self._bdf_file.write_annotation(marker)
319
320     def close_dev(self):
321         """
322             Close device connection and release resources.
323             """
324             if self.__status != iFocus.Dev.TERMINATE:
325                 # ensure socket is closed correctly
326                 self.__status = iFocus.Dev.TERMINATE_START

```

```

327         while self.__status != iFocus.Dev.TERMINATE:
328             time.sleep(0.1)
329             if self.is_alive():
330                 self.join()
331
332     def __recv_data(self):
333         try:
334             self.dev.start_data()
335             self.__status = iFocus.Dev.SIGNAL
336         except Exception:
337             self.__socket_flag = "SIGNAL mode initialization failed."
338             self.__status = iFocus.Dev.TERMINATE_START
339
340         while self.__status in [iFocus.Dev.SIGNAL]:
341             try:
342                 data = self.dev.recv_socket()
343                 if not data:
344                     raise Exception("Data transmission timeout.")
345                 ret = self.__parser.parse_data(data)
346                 if ret:
347                     if self.__lsl_eeg_flag:
348                         self._lsl_eeg.push_chunk(
349                             [frame for frame in ret for frame in frames[:-1]])
350                     )
351                     if self.__lsl_imu_flag:
352                         self._lsl_imu.push_chunk([frame[-1] for frame in ret])
353                     if self.__with_q:
354                         self.__save_data.put(ret)
355                     if self.__bdf_flag:
356                         self._bdf_file.write_chunk(ret)
357             except Exception as e:
358                 print(e)
359                 self.__socket_flag = "Data transmission timeout."
360                 self.__status = iFocus.Dev.TERMINATE_START
361
362             # clear buffer
363             self.close_lsl_eeg()
364             self.close_lsl_imu()
365             self.close_bdf_file()
366             # self.dev.stop_recv()
367             self.__parser.clear_buffer()
368             self.__save_data.put(None)
369             while self.__save_data.get() is not None:
370                 continue
371             # stop recv data
372             if self.__status != iFocus.Dev.TERMINATE_START:
373                 try: # stop data acquisition when thread ended
374                     self.dev.stop_recv()
375                 except Exception:
376                     if self.__status == iFocus.Dev.IDLE_START:
377                         self.__socket_flag = "Connection lost."
378                         self.__status = iFocus.Dev.TERMINATE_START
379
380     def run(self):
381         while self.__status != iFocus.Dev.TERMINATE_START:
382             if self.__status == iFocus.Dev.SIGNAL_START:
383                 self.__recv_data()
384             elif self.__status == iFocus.Dev.IDLE_START:
385                 self.__status = iFocus.Dev.IDLE
386                 while self.__status == iFocus.Dev.IDLE:
387                     time.sleep(0.1)
388             else:
389                 self.__socket_flag = f"Unknown status: {self.__status.name}"

```

```
        break
    try:
        self.dev.close_socket()
    finally:
        self.__status = iFocus.Dev.TERMINATE

    def __check_dev_status(self):
        if self.__socket_flag is None:
            return
        if self.is_alive():
            self.close_dev()
        raise Exception(str(self.__socket_flag))
```

3.1.1 `__init__(port=None)`

Parameters:

Name	Type	Description	Default
port	Optional[str]	if not given, connect to the first available device.	None

Source code in `src\iFocus\data_reader.py`

```

31 def __init__(self, port: Optional[str] = None) -> None:
32     """
33     Args:
34         port: if not given, connect to the first available device.
35     """
36     super().__init__(daemon=True)
37     self.__status = iFocus.Dev.TERMINATE
38     if port is None:
39         port = iFocus.find_devs()[0]
40     self.__save_data = Queue()
41     self.__parser = Parser()
42     self.dev_args = deepcopy(iFocus.dev_args)
43     self.dev = sock(port)
44     self.set_frequency()
45     self.__with_q = True
46     self.__socket_flag = "Device not connected, please connect first."
47     self.__bdf_flag = False
48     try:
49         self.dev.connect_socket()
50     except Exception as e:
51         try:
52             self.dev.close_socket()
53         finally:
54             raise e
55     self.__status = iFocus.Dev.IDLE_START
56     self.__socket_flag = None
57     self.__lsl_eeg = None
58     self.__lsl_imu = None
59     self.__lsl_imu_flag = False
60     self.__lsl_eeg_flag = False
61     self.__bdf_file = None
62     self.__enable_imu = False
63     self.dev_args["name"] = port
64     self.start()

```

3.1.2 `close_bdf_file()`

Close and save BDF file manually, invoked automatically after `stop_acquisition()` or `close_dev()`

Source code in `src\iFocus\data_reader.py`

```

284 def close_bdf_file(self):
285     """
286         Close and save BDF file manually, invoked automatically after `stop_acquisition()` or
287         `close_dev()`.
288     """
289     self.__bdf_flag = False
290     if self.__bdf_file is not None:
291         self.__bdf_file.close_bdf()
292         del self.__bdf_file

```

3.1.3 close_dev()

Close device connection and release resources.

```
Source code in src\iFocus\data_reader.py

303     def close_dev(self):
304         """
305             Close device connection and release resources.
306         """
307         if self.__status != iFocus.Dev.TERMINATE:
308             # ensure socket is closed correctly
309             self.__status = iFocus.Dev.TERMINATE_START
310             while self.__status != iFocus.Dev.TERMINATE:
311                 time.sleep(0.1)
312             if self.is_alive():
313                 self.join()
```

3.1.4 close_lsl_eeg()

Close LSL EEG stream manually, invoked automatically after `stop_acquisition()` and `close_dev()`

```
Source code in src\iFocus\data_reader.py

202     def close_lsl_eeg(self):
203         """
204             Close LSL EEG stream manually, invoked automatically after `stop_acquisition()` and
205             `close_dev()`
206         """
207         self.__lsl_eeg_flag = False
208         if hasattr(self, "_lsl_eeg"):
209             del self._lsl_eeg
```

3.1.5 close_lsl_imu()

Close LSL IMU stream manually, invoked automatically after `stop_acquisition()` and `close_dev()`

```
Source code in src\iFocus\data_reader.py

235     def close_lsl_imu(self):
236         """
237             Close LSL IMU stream manually, invoked automatically after `stop_acquisition()` and
238             `close_dev()`
239         """
240         self.__lsl_imu_flag = False
241         if hasattr(self, "_lsl_imu"):
242             del self._lsl_imu
```

3.1.6 `create_bdf_file(filename)`

Create a BDF file and save data to it, invoke it after `start_acquisition_data()`.

Parameters:

Name	Type	Description	Default
<code>filename</code>	<code>str</code>	file name to save data, accept absolute or relative path.	<code>required</code>

Raises:

Type	Description
<code>Exception</code>	if data acquisition not started or <code>save_bdf_file</code> is invoked and BDF file already created.
<code>OSError</code>	if BDF file creation failed, this may be caused by invalid file path or permission issue.
<code>importError</code>	if <code>pyedflib</code> is not installed.

Source code in `src\iFocus\data_reader.py`

```

246     def create_bdf_file(self, filename: str):
247         """
248             Create a BDF file and save data to it, invoke it after `start_acquisition_data()`.
249
250         Args:
251             filename: file name to save data, accept absolute or relative path.
252
253         Raises:
254             Exception: if data acquisition not started or `save_bdf_file` is invoked and BDF file
255             already created.
256             OSError: if BDF file creation failed, this may be caused by invalid file path or
257             permission issue.
258             ImportError: if `pyedflib` is not installed.
259         """
260         if self.__status != iFocus.Dev.SIGNAL:
261             raise Exception("Data acquisition not started")
262         if self._bdf_file is not None:
263             raise Exception("BDF file already created.")
264         from ..utils.bdfWrapper import bdfSaverEEG, bdfSaverEEGIMU
265
266         if filename[-4:].lower() != ".bdf":
267             filename += ".bdf"
268         if self.__enable_imu:
269             self._bdf_file = bdfSaverEEGIMU(
270                 filename,
271                 self.dev_args["channel_eeg"],
272                 self.dev_args["fs_eeg"],
273                 self.dev_args["channel_imu"],
274                 self.dev_args["fs_imu"],
275                 self.dev_args["type"],
276             )
277         else:
278             self._bdf_file = bdfSaverEEG(
279                 filename,
280                 self.dev_args["channel_eeg"],
281                 self.dev_args["fs_eeg"],
282                 self.dev_args["type"],
283             )
284         self.__bdf_flag = True

```

3.1.7 `find_devs()` staticmethod

Find available iFocus devices.

Returns:

Type	Description
<code>list</code>	available device ports.

Raises:

Type	Description
Exception	if no iFocus device found.

<code>Source code in src\iFocus\data_reader.py</code>	
<pre> 105 @staticmethod 106 def find_devs() -> list: 107 """ 108 Find available iFocus devices. 109 110 Returns: 111 available device ports. 112 113 Raises: 114 Exception: if no iFocus device found. 115 116 return sock._find_devs() </pre>	

3.1.8 `get_data(timeout=0.02)`

Acquire all available data, make sure this function is called in a loop when `with_q` is set to `True` in `start_acquisition_data()`

Parameters:

Name	Type	Description	Default
<code>timeout</code>	<code>Optional[float]</code>	Non-negative value, blocks at most 'timeout' seconds and return, if set to <code>None</code> , blocks until new data available.	<code>0.02</code>

Returns:

Type	Description
<code>Optional[list[Optional[list]]]</code>	A list of frames, each frame is made up of 5 eeg data and 1 imu data in a shape as below: [[eeg_0], [eeg_1], [eeg_2], [eeg_3], [eeg_4], [imu_x, imu_y, imu_z]], in which number 0~4 after _ indicates the time order of channel data.

Raises:

Type	Description
Exception	if device not connected, connection failed, data transmission timeout/init failed, or unknown error.

Data Unit

- eeg: μV
- imu: degree($^{\circ}$)

Source code in `src\eConEXG\iFocus\data_reader.py`

```

118 def get_data(
119     self, timeout: Optional[float] = 0.02
120 ) -> Optional[list[Optional[list]]]:
121     """
122         Acquire all available data, make sure this function is called in a loop when `with_q` is
123         set to `True` in `start_acquisition_data()`
124
125         Args:
126             timeout: Non-negative value, blocks at most 'timeout' seconds and return, if set to
127             `None`, blocks until new data available.
128
129         Returns:
130             A list of frames, each frame is made up of 5 eeg data and 1 imu data in a shape as
131             below:
132                 [[`eeg_0`], [`eeg_1`], [`eeg_2`], [`eeg_3`], [`eeg_4`], [`imu_x`, `imu_y`,
133                 `imu_z`]],
134                 in which number `0~4` after `_` indicates the time order of channel data.
135
136         Raises:
137             Exception: if device not connected, connection failed, data transmission timeout/init
138             failed, or unknown error.
139
140         Data Unit:
141             - eeg:  $\mu\text{V}$ 
142             - imu: degree( $^{\circ}$ )
143         """
144         self.__check_dev_status()
145         if not self.__with_q:
146             return
147         try:
148             data: list = self.__save_data.get(timeout=timeout)
149         except queue.Empty:
150             return []
151         while not self.__save_data.empty():
152             data.extend(self.__save_data.get())
153         return data

```

3.1.9 `get_dev_info()`

Get current device information, including device name, hardware channel number, acquired channels, sample frequency, etc.

Returns:

Type	Description
dict	A dictionary containing device information, which includes: <code>type</code> : hardware type; <code>channel_eeg</code> : channel dictionary, including EEG channel index and name; <code>channel_imu</code> : channel dictionary, including IMU channel index and name; <code>AdapterInfo</code> : adapter used for connection; <code>fs_eeg</code> : sample frequency of EEG in Hz; <code>fs_imu</code> : sample frequency of IMU in Hz;

Source code in `src\ConEXG\iFocus\data_reader.py`

```

90 def get_dev_info(self) -> dict:
91     """
92         Get current device information, including device name, hardware channel number, acquired
93         channels, sample frequency, etc.
94
95     Returns:
96         A dictionary containing device information, which includes:
97             `type`: hardware type;
98             `channel_eeg`: channel dictionary, including EEG channel index and name;
99             `channel_imu`: channel dictionary, including IMU channel index and name;
100            `AdapterInfo`: adapter used for connection;
101            `fs_eeg`: sample frequency of EEG in Hz;
102            `fs_imu`: sample frequency of IMU in Hz;
103        """
104     return deepcopy(self.dev_args)

```

3.1.10 `open_lsl_eeg()`

Open LSL EEG stream, can be invoked after `start_acquisition_data()`.

Raises:

Type	Description
Exception	if data acquisition not started or LSL stream already opened.
LSLException	if LSL stream creation failed.
importError	if <code>pylsl</code> is not installed or <code>liblsl</code> not installed for unix like system.

Source code in [src\iFocus\data_reader.py](#)

```

178 def open_lsl_eeg(self):
179     """
180     Open LSL EEG stream, can be invoked after `start_acquisition_data()`.
181
182     Raises:
183         Exception: if data acquisition not started or LSL stream already opened.
184         LSLError: if LSL stream creation failed.
185         ImportError: if `pylsl` is not installed or liblsl not installed for unix like
186     system.
187     """
188     if self.__status != iFocus.Dev.SIGNAL:
189         raise Exception("Data acquisition not started, please start first.")
190     if hasattr(self, "_lsl_eeg"):
191         raise Exception("LSL stream already opened.")
192     from ..utils.lslWrapper import lslSender
193
194     self._lsl_eeg = lslSender(
195         self.dev_args["channel_eeg"],
196         f"{self.dev_args['type']}EEG{self.dev_args['name'][ -2:]}",
197         "EEG",
198         self.dev_args["fs_eeg"],
199         with_trigger=False,
200     )
201     self.__lsl_eeg_flag = True

```

3.1.11 `open_lsl_imu()`

Open LSL IMU stream, can be invoked after `start_acquisition_data()`.

Raises:

Type	Description
<code>Exception</code>	if data acquisition not started or LSL stream already opened.
<code>LSLError</code>	if LSL stream creation failed.
<code>ImportError</code>	if <code>pylsl</code> is not installed or <code>liblsl</code> not installed for unix like system.

Source code in `src\iFocus\data_reader.py`

```

210 def open_lsl_imu(self):
211     """
212         Open LSL IMU stream, can be invoked after `start_acquisition_data()`.
213
214     Raises:
215         Exception: if data acquisition not started or LSL stream already opened.
216         LSLError: if LSL stream creation failed.
217         ImportError: if `pylsl` is not installed or liblsl not installed for unix like
218         system.
219     """
220     if self.__status != iFocus.Dev.SIGNAL:
221         raise Exception("Data acquisition not started, please start first.")
222     if hasattr(self, "_lsl_imu"):
223         raise Exception("LSL stream already opened.")
224     from ..utils.lslWrapper import lslSender
225
226     self._lsl_imu = lslSender(
227         self.dev_args["channel_imu"],
228         f"{self.dev_args['type']}IMU{self.dev_args['name'][ -2:]}",
229         "IMU",
230         self.dev_args["fs_imu"],
231         unit="degree",
232         with_trigger=False,
233     )
234     self.__lsl_imu_flag = True

```

3.1.12 `send_bdf_marker(marker)`

Send marker to BDF file, can be invoked after `create_bdf_file()`, otherwise it will be ignored.

Parameters:

Name	Type	Description	Default
<code>marker</code>	<code>str</code>	marker string to write.	<i>required</i>

Source code in `src\iFocus\data_reader.py`

```

293 def send_bdf_marker(self, marker: str):
294     """
295         Send marker to BDF file, can be invoked after `create_bdf_file()`, otherwise it will be
296         ignored.
297
298     Args:
299         marker: marker string to write.
300     """
301     if hasattr(self, "_bdf_file"):
302         self._bdf_file.write_annotation(marker)

```

3.1.13 set_frequency(fs_eeg=None)

Change the sampling frequency of iFocus.

Parameters:

Name	Type	Description	Default
fs_eeg	int	sampling frequency of eeg data, should be 250 or 500, fs_imu will be automatically set to 1/5 of fs_eeg.	None

Raises:

Type	Description
ValueError	if fs_eeg is not 250 or 500.
NotImplementedError	device firmware too old, not supporting 500Hz.

Source code in [src\ConEXG\iFocus\data_reader.py](#)

```

66 def set_frequency(self, fs_eeg: int = None):
67     """
68     Change the sampling frequency of iFocus.
69
70     Args:
71         fs_eeg: sampling frequency of eeg data, should be 250 or 500,
72                 fs_imu will be automatically set to 1/5 of fs_eeg.
73
74     Raises:
75         ValueError: if fs_eeg is not 250 or 500.
76         NotImplementedError: device firmware too old, not supporting 500Hz.
77     """
78     if self.__status == iFocus.Dev.SIGNAL:
79         raise Exception("Data acquisition already started, please stop first.")
80     if fs_eeg is None:
81         fs_eeg = self.dev_args["fs_eeg"]
82     if fs_eeg not in [250, 500]:
83         raise ValueError("fs_eeg should be 250 or 500")
84     self.dev_args["fs_eeg"] = fs_eeg
85     fs_imu = fs_eeg // 5
86     self.dev_args["fs_imu"] = fs_imu
87     if hasattr(self, "dev"):
88         self.dev.set_frequency(fs_eeg)

```

3.1.14 start_acquisition_data(with_q=True)

Send data acquisition command to device, block until data acquisition started or failed.

Parameters:

Name	Type	Description	Default
with_q	bool	if True, signal data will be stored in a queue and should be acquired by calling <code>get_data()</code> in a loop in case data queue is full. if False, new data will not be directly available and can only be acquired through lsI stream.	True

Source code in `src\iFocus\data_reader.py`

```

150 def start_acquisition_data(self, with_q: bool = True) -> None:
151     """
152         Send data acquisition command to device, block until data acquisition started or failed.
153
154     Args:
155         with_q: if True, signal data will be stored in a queue and **should** be acquired by
156         calling `get_data()` in a loop in case data queue is full.
157             if False, new data will not be directly available and can only be acquired
158         through lsI stream.
159
160     """
161     self.__check_dev_status()
162     self.__with_q = with_q
163     if self.__status == iFocus.Dev.SIGNAL:
164         return
165     self.__status = iFocus.Dev.SIGNAL_START
166     while self.__status not in [iFocus.Dev.SIGNAL, iFocus.Dev.TERMINATE]:
167         time.sleep(0.01)
168     self.__check_dev_status()

```

3.1.15 `stop_acquisition()`

Stop data or impedance acquisition, block until data acquisition stopped or failed.

Source code in `src\iFocus\data_reader.py`

```

168 def stop_acquisition(self) -> None:
169     """
170         Stop data or impedance acquisition, block until data acquisition stopped or failed.
171     """
172     self.__check_dev_status()
173     self.__status = iFocus.Dev.IDLE_START
174     while self.__status not in [iFocus.Dev.IDLE, iFocus.Dev.TERMINATE]:
175         time.sleep(0.01)
176     self.__check_dev_status()

```

I.I Trigger Box

4 Wireless

Source code in `src\triggerBox\triggerbox.py`

```

5  class triggerBoxWireless:
6      def __init__(self, port: str = None):
7          """
8              Args:
9                  port: The serial port of the trigger box. If not given,
10                     the function will try to find the trigger box automatically.
11
12             Raises:
13                 Exception: If the trigger box is not found.
14
15         from serial import Serial
16         from serial.tools.list_ports import comports
17
18     if not port:
19         for ports in comports():
20             if ports.pid == 0x6001 and ports.vid == 0x0403:
21                 port = ports.device
22                 break
23             else:
24                 raise Exception("Trigger box not found")
25         self.dev = Serial(port, baudrate=115200, timeout=1)
26         self.__last_timestamp = time.perf_counter()
27         self.__warn = "Marker interval too short, amplifier may fail to receive it. Suggested
28 interval is above 50ms"
29         time.sleep(0.1)
30
31     def sendMarker(self, marker: int):
32         """
33             Send a marker to the trigger box.
34
35             Args:
36                 marker: range from `1` to `255`, `13` is not available and reserved for internal
37 use.
38
39             Raises:
40                 Exception: If the marker is invalid.
41
42             if time.perf_counter() - self.__last_timestamp < 0.04:
43                 print(self.__warn)
44             if not isinstance(marker, int):
45                 marker = int(marker)
46             if marker == 13 or marker <= 0 or marker > 255:
47                 raise Exception("Invalid marker")
48             marker = marker.to_bytes(length=1, byteorder="big", signed=False)
49             self.dev.write(marker + b"\x55\x66\x0d")
50             self.__last_timestamp = time.perf_counter()
51
52     def close_dev(self):
53         self.dev.close()

```

4.1 __init__(port=None)

Parameters:

Name	Type	Description	Default
port	str	The serial port of the trigger box. If not given, the function will try to find the trigger box automatically.	None

Raises:

Type	Description
Exception	If the trigger box is not found.

Source code in [src\triggerBox\triggerbox.py](#) ▾

```

6 def __init__(self, port: str = None):
7     """
8     Args:
9         port: The serial port of the trigger box. If not given,
10            the function will try to find the trigger box automatically.
11
12    Raises:
13        Exception: If the trigger box is not found.
14    """
15    from serial import Serial
16    from serial.tools.list_ports import comports
17
18    if not port:
19        for ports in comports():
20            if ports.pid == 0x6001 and ports.vid == 0x0403:
21                port = ports.device
22                break
23        else:
24            raise Exception("Trigger box not found")
25    self.dev = Serial(port, baudrate=115200, timeout=1)
26    self.__last_timestamp = time.perf_counter()
27    self.__warn = "Marker interval too short, amplifier may fail to receive it. Suggested
28    interval is above 50ms"
29    time.sleep(0.1)

```

4.2 sendMarker(marker)

Send a marker to the trigger box.

Parameters:

Name	Type	Description	Default
marker	int	range from 1 to 255, 13 is not available and reserved for internal use.	required

Raises:

Type	Description
Exception	If the marker is invalid.

Source code in [src\ConEXG\triggerBox\triggerbox.py](#) ▾

```

30 def sendMarker(self, marker: int):
31     """
32     Send a marker to the trigger box.
33
34     Args:
35         marker: range from `1` to `255`, `13` is not available and reserved for internal use.
36
37     Raises:
38         Exception: If the marker is invalid.
39     """
40     if time.perf_counter() - self.__last_timestamp < 0.04:
41         print(self.__warn)
42     if not isinstance(marker, int):
43         marker = int(marker)
44     if marker == 13 or marker <= 0 or marker > 255:
45         raise Exception("Invalid marker")
46     marker = marker.to_bytes(length=1, byteorder="big", signed=False)
47     self.dev.write(marker + b"\x55\x66\x0d")
48     self.__last_timestamp = time.perf_counter()

```

5 Wired

Source code in `src\triggerBox\triggerbox.py`

```
54 class triggerBoxWired:
55     def __init__(self, port: str = None):
56         """
57             Args:
58                 port: The serial port of the trigger box. If not given,
59                     the function will try to find the trigger box automatically.
60
61             Raises:
62                 Exception: If the trigger box is not found.
63         """
64         from serial import Serial
65         from serial.tools.list_ports import comports
66
67         if not port:
68             for ports in comports():
69                 if ports.pid == 0x5740 and ports.vid == 0x0483:
70                     port = ports.device
71                     break
72             else:
73                 raise Exception("Trigger box not found")
74         self.dev = Serial(port, timeout=1)
75
76     def sendMarker(self, marker: int):
77         """
78             Send a marker to the trigger box.
79
80             Args:
81                 marker: range from `1` to `255`.
82
83             Raises:
84                 Exception: If the marker is invalid.
85         """
86         if not isinstance(marker, int):
87             marker = int(marker)
88         if marker <= 0 or marker > 255:
89             raise Exception("Invalid marker")
90         self.dev.write(marker.to_bytes(length=1, byteorder="big", signed=False))
91
92     def close_dev(self):
93         self.dev.close()
```

5.1 __init__(port=None)

Parameters:

Name	Type	Description	Default
port	str	The serial port of the trigger box. If not given, the function will try to find the trigger box automatically.	None

Raises:

Type	Description
Exception	If the trigger box is not found.

Source code in [src\triggerBox\triggerbox.py](#)

```

55 def __init__(self, port: str = None):
56     """
57     Args:
58         port: The serial port of the trigger box. If not given,
59             the function will try to find the trigger box automatically.
60
61     Raises:
62         Exception: If the trigger box is not found.
63     """
64     from serial import Serial
65     from serial.tools.list_ports import comports
66
67     if not port:
68         for ports in comports():
69             if ports.pid == 0x5740 and ports.vid == 0x0483:
70                 port = ports.device
71                 break
72             else:
73                 raise Exception("Trigger box not found")
74     self.dev = Serial(port, timeout=1)

```

5.2 sendMarker(marker)

Send a marker to the trigger box.

Parameters:

Name	Type	Description	Default
marker	int	range from 1 to 255.	required

Raises:

Type	Description
Exception	If the marker is invalid.

Source code in [src\eConEXG\triggerBox\triggerbox.py](#) ▾

```
76 def sendMarker(self, marker: int):
77     """
78     Send a marker to the trigger box.
79
80     Args:
81         marker: range from `1` to `255`.
82
83     Raises:
84         Exception: If the marker is invalid.
85     """
86     if not isinstance(marker, int):
87         marker = int(marker)
88     if marker <= 0 or marker > 255:
89         raise Exception("Invalid marker")
90     self.dev.write(marker.to_bytes(length=1, byteorder="big", signed=False))
```

6 Light Stimulator

Source code in [src\eConEXG\triggerBox\triggerbox.py](#)

```

96  class lightStimulator:
97      def __init__(self, port: str = None):
98          from serial import Serial
99          from serial.tools.list_ports import comports
100
101         self.wait_time = 0.1
102         self.channels = 6
103
104         if not port:
105             for ports in comports():
106                 if (
107                     ports.pid == 0x6001
108                     and ports.vid == 0x0403
109                     and ports.serial_number in ["LIGHTSTIMA", "LIGHTSTIM"]
110                 ):
111                     port = ports.device
112                     break
113                 else:
114                     raise Exception("Light stimulator not found")
115         self.dev = Serial(port, baudrate=115200, timeout=2)
116         self.dev.read_all()
117
118     def vep_mode(self, fs: list[Optional[float]] = [1, 1, 1, 1, 1, 1]):
119         """
120             Enter VEP mode, which allows you to control the frequency of each channel separately.
121
122             Args:
123                 fs: List of frequencies in Hz, range from 0 to 100 with 0.1Hz resolution.
124                 If a corresponding frequency is None, 0 or not given, it will be set to off.
125
126             Raises:
127                 Exception: If the frequency is invalid or hardware error.
128         """
129         fss = fs.copy()
130         if len(fss) < self.channels:
131             fss += [0] * (self.channels - len(fss))
132         for i in range(self.channels):
133             fss[i] = self._validate_fs(fss[i])
134         command = ",".join([f"{f:.1f}" for f in fss[: self.channels]])
135         command = f"AT+VEP={command}\r\n".encode()
136         self.dev.write(command)
137         time.sleep(self.wait_time)
138         ret = self.dev.read_all()
139         if b"SSVEP MODE OK" not in ret:
140             raise Exception("Failed to set VEP mode")
141
142     def erp_mode(self, fs: float):
143         """
144             Enter ERP mode, which allows you to control the frequency of all channels at once.
145
146             Args:
147                 fs: Frequency in Hz, range from 0 to 100 with 0.1Hz resolution.
148
149             Raises:
150                 Exception: If the frequency is invalid or hardware error.
151         """
152         fs = self._validate_fs(fs)
153         command = f"AT+ERP={fs:.1f}\r\n".encode()
154         self.dev.write(command)
155         time.sleep(self.wait_time)
156         ret = self.dev.read_all()
157         if b"ERP MODE OK" not in ret:

```

```

158         raise Exception("Failed to set VEP mode")
159
160     def _validate_fs(self, fs: Optional[float]):
161         if not isinstance(fs, (int, float)):
162             if fs is None:
163                 fs = 0
164             else:
165                 raise Exception("Invalid frequency")
166         if fs > 100 or fs < 0:
167             raise Exception("Invalid frequency")
168         return fs
169
170     def close_dev(self):
171         self.dev.close()

```

6.1 erp_mode(fs)

Enter ERP mode, which allows you to control the frequency of all channels at once.

Parameters:

Name	Type	Description	Default
fs	float	Frequency in Hz, range from 0 to 100 with 0.1Hz resolution.	<i>required</i>

Raises:

Type	Description
Exception	If the frequency is invalid or hardware error.

Source code in `src\triggerBox\triggerbox.py`

```

142     def erp_mode(self, fs: float):
143         """
144             Enter ERP mode, which allows you to control the frequency of all channels at once.
145
146             Args:
147                 fs: Frequency in Hz, range from 0 to 100 with 0.1Hz resolution.
148
149             Raises:
150                 Exception: If the frequency is invalid or hardware error.
151             """
152             fs = self._validate_fs(fs)
153             command = f"AT+ERP={fs:.1f}\r\n".encode()
154             self.dev.write(command)
155             time.sleep(self.wait_time)
156             ret = self.dev.read_all()
157             if b"ERP MODE OK" not in ret:
158                 raise Exception("Failed to set VEP mode")

```

6.2 vep_mode(fs=[1, 1, 1, 1, 1, 1])

Enter VEP mode, which allows you to control the frequency of each channel separately.

Parameters:

Name	Type	Description	Default
fs	list[Optional[float]]	List of frequencies in Hz, range from 0 to 100 with 0.1Hz resolution. If a corresponding frequency is None, 0 or not given, it will be set to off.	[1, 1, 1, 1, 1, 1]

Raises:

Type	Description
Exception	If the frequency is invalid or hardware error.

Source code in [src\eConEXG\triggerBox\triggerbox.py](#)

```
118 def vep_mode(self, fs: list[Optional[float]] = [1, 1, 1, 1, 1, 1]):
119     """
120     Enter VEP mode, which allows you to control the frequency of each channel separately.
121
122     Args:
123         fs: List of frequencies in Hz, range from 0 to 100 with 0.1Hz resolution.
124             If a corresponding frequency is None, 0 or not given, it will be set to off.
125
126     Raises:
127         Exception: If the frequency is invalid or hardware error.
128     """
129     fss = fs.copy()
130     if len(fss) < self.channels:
131         fss += [0] * (self.channels - len(fss))
132     for i in range(self.channels):
133         fss[i] = self._validate_fs(fss[i])
134     command = ",".join([f"{f:.1f}" for f in fss[: self.channels]])
135     command = f"AT+VEP={command}\r\n".encode()
136     self.dev.write(command)
137     time.sleep(self.wait_time)
138     ret = self.dev.read_all()
139     if b"SSVEP MODE OK" not in ret:
140         raise Exception("Failed to set VEP mode")
```

7 Changelog

Here you can find all the released changes to eConEXG.

7.1 0.1.18

Release on 2024-10-10

- **ADD** Add `create_bdf_file` and `close_bdf_file` to `iRecorder` and `eConAlpha`.
- **ADD** Add imu support to `iRecorder` and `eConAlpha`.

7.2 0.1.15

- **ADD** `eConAlpha` device SDK

7.3 0.1.14

Released on 2024-08-08.

- **Update** `pyEDFlib` dependency requirement in `pyproject.toml` from `0.1.37` to `0.1.38` to support `numpy>=2.0.0`.
- **Optimize** `set_frequency()` in `iRecorder`, now you can set sample rate after device connection.
- **Fix** `iFocus` not raise Exception after lost connection with `USBAdapter`.
- **Change** `with_q` argument in `iRecorder` and `iFocus` from constructor to `start_acquisition_data()`.
- **Deprecate** `save_bdf_file()` in `iRecorder`, use `create_bdf_file()` instead.

7.4 0.1.13

Released on 2024-08-01.

- **Add** `get_dev_info()` in `iFocus` class to get device information.
- **Add** selectable `500Hz` `eeg` and corresponding `100Hz` IMU sampling rate in `iFocus` class through `set_frequency()`.
- **Add** `with_q = False` option in `iRecorder` and `iFocus` constructor to drop the necessity of loop calling `get_data()` in SIGNAL mode.
- **Add** `__version__` field of `eConEXG` package to check the package version, it can be accessed through `eConEXG.__version__`.

- **Fix** last valid packet number in `iFocus` warning message wrongly displayed as a fixed number issue.
- **Change** the default data parse length from 0.01 seconds to 10 frames in `iRecorder` to match hardware settings.
- **Improve** the aesthetics of a document interface.

7.5 0.1.12

Released on 2024-07-25.

- **Fix** equipment format warning issue on bdf save.
- **Fix** Rounding `physical_max` value of bdf file, resulting to more accurate data precision.
- **Fix** `sendMarker()` function not working in `triggerBox` class when python version<3.11.

7.6 0.1.11

Released on 2024-07-12.

- **Update** documentation homepage.

7.7 0.1.10

Released on 2024-07-08.

- **Add** `Isl` support for `iFocus`.
 - **Add** support for `iRecorder` 16-channel wired mode.
-
- **Change** default timeout of `get_data()` function in `iRecorder` from `None` to `0.02`.

7.8 0.1.09

Released on 2024-06-28.

- **Add** support for `iRecorder` 8-channel wired mode.